

MONIDOR



Monidor Vitals Remote Patient Monitoring

For continuous monitoring of SpO₂ and pulse rate in hospital wards.



Detect Issues Before They Escalate

Monidor Vitals™ alerts if the patient's SpO₂ or pulse rate deviates from the limits, as well as for technical reasons. Alarm limits can be set per patient.



"Easy to use and quick to start monitoring."

"The overall picture of oxygenation is better understood. Draws attention if changes occur. Handy with isolation patients."

**Internal medicine ward 4A,
Savonlinna Central Hospital**

Operating Principle

- 1 The pulse oximeter continuously measures the patient's oxygen saturation and pulse.
- 2 The pulse oximeter is connected via Bluetooth to a base station in the patient's room, which transmits the readings via Wi-Fi to a remote monitoring device.
- 3 Monitor the patient's oxygen saturation and pulse remotely using a PC or smartphone.

Monidor Vitalis remote monitoring application and Nonin WristOx2 pulse oximeter are CE-marked medical devices.

MONIDOR